

NEEP21 Polyclonal Antibody

Cat No: HR1AP4025

For research use only

Overview

Product Name	NEEP21 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	NSG1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NSG1; D4S234; Neuron-specific protein family member 1; Brain neuron cytoplasmic protein 1
Background	similarity:Belongs to the NSG family.,