

Cables1 Polyclonal Antibody

Cat No: HR1AP2931

For research use only

Overview

Product Name	Cables1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	67kDa
GeneID?Human?	CABLES1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CABLES1; CABLES; CDK5 and ABL1 enzyme substrate 1; Interactor with CDK3 1; Ik3-1
Background	Cdk5 and Abl enzyme substrate 1(CABLES1) Homo sapiens This gene encodes a protein involved in regulation of the cell cycle through interactions with several cyclin-dependent kinases. One study (PMID: 16177568) reported aberrant splicing of transcripts from this gene which results in removal of the cyclin binding domain only in human cancer cells, and reduction in gene expression was shown in colorectal cancers (PMID: 17982127).Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2012],