

DOR-1 (phospho Ser363) Polyclonal Antibody

Cat No: HR1AP7552

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | DOR-1 (phospho Ser363) Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IHC-p,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 36kDa |
| GeneID?Human? | OPRD1 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | OPRD1; OPRD; Delta-type opioid receptor; D-OR-1; DOR-1 |
| Background | function:Inhibits neurotransmitter release by reducing calcium ion currents and increasing potassium ion conductance. Highly stereoselective. receptor for enkephalins.,online information:Delta opioid receptor entry,similarity:Belongs to the G-protein coupled receptor 1 family.,subunit:Interacts with GPRASP1., |