

GRIN1 Polyclonal Antibody

Cat No: HR1AP3585

For research use only

Overview

Product Name	GRIN1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	102kDa
GeneID?Human?	GPRIN1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPRIN1; KIAA1893; G protein-regulated inducer of neurite outgrowth 1; GRIN1
Background	function:May be involved in neurite outgrowth.,PTM:Palmitoylation on Cys-999 and/or Cys-1000 is required for membrane targeting.,subcellular location:Highly enriched in growth cone.,subunit:Interacts with activated forms of GNAI1, GNAO1 and GNAZ.,tissue specificity:Widely expressed in the central nervous system, with highest levels in spinal cord.,