

GABA T-1 Polyclonal Antibody

Cat No: HR1AP9698

For research use only

Overview

Product Name	GABA T-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	SLC6A1 GABATR GABT1 GAT1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Sodium- and chloride-dependent GABA transporter 1 (GAT-1) (Solute carrier family 6 member 1)
Background	solute carrier family 6 member 1(SLC6A1) Homo sapiens The SLC6A1 gene encodes a gamma-aminobutyric acid (GABA) transporter, which removes GABA from the synaptic cleft (Hirunsatit et al., 2009 [PubMed 19077666]).[supplied by OMIM, Jul 2009],