

DOR-1 Polyclonal Antibody

Cat No: HR1AP7554

For research use only

Overview

Product Name	DOR-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	OPRD1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	OPRD1; OPRD; Delta-type opioid receptor; D-OR-1; DOR-1
Background	function:Inhibits neurotransmitter release by reducing calcium ion currents and increasing potassium ion conductance. Highly stereoselective. receptor for enkephalins.,online information:Delta opioid receptor entry,similarity:Belongs to the G-protein coupled receptor 1 family.,subunit:Interacts with GPRASP1.,