

## Stra8 Polyclonal Antibody

Cat No: HR1AP5398

For research use only

### Overview

Product Name	Stra8 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	66kDa
GeneID?Human?	STRA8
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	STRA8; Stimulated by retinoic acid gene 8 protein homolog
Background	<p>stimulated by retinoic acid 8(STRA8) Homo sapiens This gene encodes a retinoic acid-responsive protein. A homologous protein in mouse has been shown to be involved in the regulation of meiotic initiation in both spermatogenesis and oogenesis, though feature differences between the mouse and human proteins suggest that these homologs are not entirely functionally equivalent. It is thought that this gene may play a role in spermatogenesis in humans. [provided by RefSeq, Nov 2010],</p>