

Synuclein-? Polyclonal Antibody

Cat No: HR1AP4652

For research use only

Overview

Product Name	Synuclein-? Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	14kDa
GeneID?Human?	SNCB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SNCB; Beta-synuclein
Background	synuclein beta(SNCB) Homo sapiens This gene encodes a member of a small family of proteins that inhibit phospholipase D2 and may function in neuronal plasticity. The encoded protein is abundant in lesions of patients with Alzheimer disease. A mutation in this gene was found in individuals with dementia with Lewy bodies. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2015],