

TMEM37 Polyclonal Antibody

Cat No: HR1AP5939

For research use only

Overview

Product Name	TMEM37 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	TMEM37
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TMEM37; PR; Voltage-dependent calcium channel gamma-like subunit; Neuronal voltage-gated calcium channel gamma-like subunit; Transmembrane protein 37
Background	function:Thought to stabilize the calcium channel in an inactivated (closed) state. Modulates calcium current when coexpressed with CACNA1G.,similarity:Belongs to the PMP-22/EMP/MP20 family. CACNG subfamily.,subunit:The L-type calcium channel is composed of five subunits: alpha-1, alpha-2/delta, beta and gamma.,