

# TMEM101 Polyclonal Antibody

Cat No: HR1AP8861

For research use only

## Overview

Product Name	TMEM101 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	TMEM101
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TMEM101; Transmembrane protein 101; Putative NF-kappa-B-activating protein 130
Background	function:May activate NF-kappa-B signaling pathways.,