

Relaxin Receptor 4 Polyclonal Antibody

Cat No: HR1AP6926

For research use only

Overview

Product Name	Relaxin Receptor 4 Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	RXFP4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RXFP4; GPR100; RLN3R2; Relaxin-3 receptor 2; RLN3 receptor 2; G-protein coupled receptor 100; G-protein coupled receptor GPCR142; Insulin-like peptide INSL5 receptor; Relaxin family peptide receptor 4
Background	relaxin/insulin like family peptide receptor 4(RXFP4) Homo sapiens GPR100 is a member of the rhodopsin family of G protein-coupled receptors (GPRs) (Fredriksson et al., 2003 [PubMed 14623098]).[supplied by OMIM, Mar 2008],