

Olfactory receptor 2J3 Polyclonal Antibody

Cat No: HR1AP7387

For research use only

Overview

Product Name	Olfactory receptor 2J3 Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	OR2J3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	OR2J3; Olfactory receptor 2J3; Hs6M1-3; Olfactory receptor OR6-16; OR6-6; Olfactory receptor 6-6
Background	olfactory receptor family 2 subfamily J member 3(OR2J3) Homo sapiens This gene encodes a G-protein-coupled receptor (GPCR) that functions as an olfactory receptor. Olfactory receptors interact with odorant molecules in the nose to initiate a neuronal response that triggers the perception of a smell. The protein encoded by this gene responds to cis-3-hexen-1-ol, which is released by wounded plants, including cut grass. This gene is situated in a cluster of similar olfactory-receptor coding genes on chromosome 6. [provided by RefSeq, May 2013],