

## Mas1-like Polyclonal Antibody

Cat No: HR1AP5772

For research use only

### Overview

Product Name	Mas1-like Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	MAS1L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MAS1L; MRG; Mas-related G-protein coupled receptor MRG; MAS-R; MAS1-like
Background	similarity:Belongs to the G-protein coupled receptor 1 family.,similarity:Belongs to the G-protein coupled receptor 1 family. Mas subfamily.,