

LDLRAD3 Polyclonal Antibody

Cat No: HR1AP5950

For research use only

Overview

Product Name	LDLRAD3 Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	LDLRAD3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	LDLRAD3; Low-density lipoprotein receptor class A domain-containing protein 3
Background	similarity:Belongs to the LDLR family.,similarity:Contains 3 LDL-receptor class A domains.,