

# LRP11 Polyclonal Antibody

Cat No: HR1AP8897

For research use only

## Overview

Product Name	LRP11 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	LRP11
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	LRP11; Low-density lipoprotein receptor-related protein 11; LRP-11
Background	similarity:Belongs to the LDLR family.,similarity:Contains 1 LDL-receptor class A domain.,similarity:Contains 1 MANSC domain.,similarity:Contains 1 PKD domain.,