

LRP10 Polyclonal Antibody

Cat No: HR1AP6570

For research use only

Overview

Product Name	LRP10 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	76kDa
GeneID?Human?	LRP10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	LRP10; MSTP087; SP220; Low-density lipoprotein receptor-related protein 10; LRP-10
Background	LDL receptor related protein 10(LRP10) Homo sapiens This gene encodes a low density lipoprotein receptor family protein. A similar protein in mouse is thought to play a role in the uptake of apolipoprotein E-containing lipoproteins. [provided by RefSeq, Jul 2016],