

GPRC5D Polyclonal Antibody

Cat No: HR1AP7840

For research use only

Overview

Product Name	GPRC5D Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	GPRC5D
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPRC5D; G-protein coupled receptor family C group 5 member D
Background	G protein-coupled receptor class C group 5 member D(GPRC5D) Homo sapiens The protein encoded by this gene is a member of the G protein-coupled receptor family; however, the specific function of this gene has not yet been determined. [provided by RefSeq, Jul 2008],