

GPR42 Polyclonal Antibody

Cat No: HR1AP6730

For research use only

Overview

Product Name	GPR42 Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	GPR42
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR42; GPR42P; G-protein coupled receptor 42
Background	GPR42, also named as GPR42P, FFAR1L and GPR41L, belongs to the G-protein coupled receptor 1 family. GPR42 appears to have arisen by tandem duplication of FFAR3 and have lost the ability to activate G-protein family in response to carboxylate ligands primarily due to an amino acid change at position 174. It is an orphan receptor.