

## GPR63 Polyclonal Antibody

Cat No: HR1AP8751

For research use only

### Overview

|                       |   |
|-----------------------|---|
| Product Name          | GPR63 Polyclonal Antibody   |
| Source                | Rabbit  |
| Applications          | IF,ELISA  |
| Species Reactivity    | Human,Mouse   |
| Recommended Dilutions |   |
| Immunogen             |   |
| Species               | Rabbit  |
| Storage               | -20°C/1 year  |
| Isotype               |   |
| Clonality             |   |
| Concentration         | 1 mg/ml   |
| Observed band         | kDa   |
| GeneID?Human?         | GPR63   |
| Human Swiss-Prot No.  |   |
| Cellular localization |   |
| Alternative Names     | GPR63; PSP24B; Probable G-protein coupled receptor 63; PSP24-2; PSP24-beta  |
| Background            | G protein-coupled receptor 63(GPR63) Homo sapiens This gene encodes a G protein-coupled receptor. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Dec 2011], |