

GPR43 Polyclonal Antibody

Cat No: HR1AP6731

For research use only

Overview

Product Name	GPR43 Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	FFAR2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FFAR2; GPR43; Free fatty acid receptor 2; G-protein coupled receptor 43
Background	free fatty acid receptor 2(FFAR2) Homo sapiens This gene encodes a member of the GP40 family of G protein-coupled receptors that are clustered together on chromosome 19. The encoded protein is a receptor for short chain free fatty acids and may be involved in the inflammatory response and in regulating lipid plasma levels. [provided by RefSeq, Apr 2009],