

GPR31 Polyclonal Antibody

Cat No: HR1AP6717

For research use only

Overview

Product Name	GPR31 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	GPR31
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR31; 12-(S)-hydroxy-5; 8,10,14-eicosatetraenoic acid receptor; 12-(S)-HETE receptor; 12-HETER; G-protein coupled receptor 31
Background	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,