

## GPR25 Polyclonal Antibody

Cat No: HR1AP6713

For research use only

### Overview

Product Name	GPR25 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	GPR25
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR25; Probable G-protein coupled receptor 25
Background	G protein-coupled receptor 25(GPR25) Homo sapiens This gene is intronless and encodes a member of the G-protein coupled receptor 1 family. G-protein coupled receptors are membrane proteins which activate signaling cascades as a response to extracellular stress. This gene has been linked to arterial stiffness. [provided by RefSeq, Nov 2012],