

GPR20 Polyclonal Antibody

Cat No: HR1AP6706

For research use only

Overview

Product Name	GPR20 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	GPR20
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR20; G-protein coupled receptor 20
Background	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,tissue specificity:Liver and the brain regions thalamus, putamen, and caudate, but not in frontal cortex, pons and hypothalamus.,