

## GPR160 Polyclonal Antibody

Cat No: HR1AP6636

For research use only

### Overview

Product Name	GPR160 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	GPR160
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR160; GPCR150; Probable G-protein coupled receptor 160; G-protein coupled receptor GPCR1; hGPCR1
Background	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,