

GPR153 Polyclonal Antibody

Cat No: HR1AP7148

For research use only

Overview

Product Name	GPR153 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	65kDa
GeneID?Human?	GPR153
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR153; PGR1; Probable G-protein coupled receptor 153; G-protein coupled receptor PGR1
Background	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,