

GPR151 Polyclonal Antibody

Cat No: HR1AP5884

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | GPR151 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IF,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 46kDa |
| GeneID?Human? | GPR151 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | GPR151; PGR7; Probable G-protein coupled receptor 151; G-protein coupled receptor PGR7; GPCR-2037 |
| Background | function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family., |