

## GPR144 Polyclonal Antibody

Cat No: HR1AP6953

For research use only

### Overview

Product Name	GPR144 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	104kDa
GeneID?Human?	GPR144
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR144; PGR24; Probable G-protein coupled receptor 144; G-protein coupled receptor PGR24
Background	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily.,similarity:Contains 1 GPS domain.,similarity:Contains 1 pentaxin domain.,