

GALR3 Polyclonal Antibody

Cat No: HR1AP8894

For research use only

Overview

Product Name	GALR3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	GALR3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GALR3; GALNR3; Galanin receptor type 3; GAL3-R; GALR-3
Background	galanin receptor 3(GALR3) Homo sapiens The neuropeptide galanin modulates a variety of physiologic processes including cognition/memory, sensory/pain processing, hormone secretion, and feeding behavior. The human galanin receptors are G protein-coupled receptors that functionally couple to their intracellular effector through distinct signaling pathways. GALR3 is found in many tissues and may be expressed as 1.4-, 2.4-, and 5-kb transcripts [provided by RefSeq, Jul 2008].