

GPR1 Polyclonal Antibody

Cat No: HR1AP6685

For research use only

Overview

Product Name	GPR1 Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	GPR1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR1; G-protein coupled receptor 1
Background	function:Orphan receptor. Probable opioid receptor. Can act as a co-receptor for HIV-1.,similarity:Belongs to the G-protein coupled receptor 1 family.,tissue specificity:Expressed in hippocampus.,