

Peropsin Polyclonal Antibody

Cat No: HR1AP5646

For research use only

Overview

Product Name	Peropsin Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	RRH
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RRH; Visual pigment-like receptor peropsin
Background	retinal pigment epithelium-derived rhodopsin homolog(RRH) Homo sapiens Opsins are members of the guanine nucleotide-binding protein (G protein)-coupled receptor superfamily. This gene belongs to the seven-exon subfamily of mammalian opsin genes that includes opsin 5 and retinal G protein coupled receptor. [provided by RefSeq, Jul 2008],