

Olfactory receptor 8J3 Polyclonal Antibody

Cat No: HR1AP4186

For research use only

Overview

Product Name	Olfactory receptor 8J3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	OR8J3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	OR8J3; Olfactory receptor 8J3; Olfactory receptor OR11-173
Background	function:Odorant receptor .,similarity:Belongs to the G-protein coupled receptor 1 family.,