

# LDLRAD2 Polyclonal Antibody

Cat No: HR1AP3818

For research use only

## Overview

Product Name	LDLRAD2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	LDLRAD2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	LDLRAD2; Low-density lipoprotein receptor class A domain-containing protein 2
Background	similarity:Belongs to the LDLR family.,similarity:Contains 1 LDL-receptor class A domain.,