

LDLRAD1 Polyclonal Antibody

Cat No: HR1AP3817

For research use only

Overview

Product Name	LDLRAD1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	LDLRAD1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	LDLRAD1; Low-density lipoprotein receptor class A domain-containing protein 1
Background	similarity:Belongs to the LDLR family.,similarity:Contains 3 LDL-receptor class A domains.,