

JTB Polyclonal Antibody

Cat No: HR1AP5679

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | JTB Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 17kDa |
| GeneID?Human? | JTB |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | JTB; HSPC222; Protein JTB; Jumping translocation breakpoint protein; Prostate androgen-regulated protein; PAR protein |
| Background | similarity:Belongs to the JTB family.,tissue specificity:Expressed in all normal human tissues studied but overexpressed in most of their malignant counterparts., |