

GPR75 Polyclonal Antibody

Cat No: HR1AP5684

For research use only

Overview

Product Name	GPR75 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	59kDa
GeneID?Human?	GPR75
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR75; Probable G-protein coupled receptor 75
Background	G protein-coupled receptor 75(GPR75) Homo sapiens GPR75 is a member of the G protein-coupled receptor family. GPRs are cell surface receptors that activate guanine-nucleotide binding proteins upon the binding of a ligand.[supplied by OMIM, Jul 2002],