

## GPR19 Polyclonal Antibody

Cat No: HR1AP3572

For research use only

### Overview

Product Name	GPR19 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GeneID?Human?	GPR19
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR19; Probable G-protein coupled receptor 19; GPR-NGA
Background	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,tissue specificity:Abundant expression in the brain.,