

GPR38 Polyclonal Antibody

Cat No: HR1AP3574

For research use only

Overview

Product Name	GPR38 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	MLNR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MLNR; GPR38; MTLR; MTLR1; Motilin receptor; G-protein coupled receptor 38
Background	<p>motilin receptor(MLNR) Homo sapiens Motilin is a 22 amino acid peptide hormone expressed throughout the gastrointestinal (GI) tract. The protein encoded by this gene is a motilin receptor which is a member of the G-protein coupled receptor 1 family. This member is a multi-pass transmembrane protein, and is an important therapeutic target for the treatment of hypomotility disorders. [provided by RefSeq, Aug 2011],</p>