

GPR114 Polyclonal Antibody

Cat No: HR1AP6336

For research use only

Overview

Product Name	GPR114 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	59kDa
GeneID?Human?	GPR114
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR114; PGR27; Probable G-protein coupled receptor 114; G-protein coupled receptor PGR27
Background	adhesion G protein-coupled receptor G5(ADGRG5) Homo sapiens This gene encodes a member of the adhesion family of G-protein coupled receptors. Members of this family are characterized by long N-termini and multiple functional domains. They may play a role in the immune system as well as in the central nervous system. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016],