

GPR107 Polyclonal Antibody

Cat No: HR1AP3557

For research use only

Overview

Product Name	GPR107 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	67kDa
GeneID?Human?	GPR107
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR107; KIAA1624; LUSTR1; Protein GPR107; Lung seven transmembrane receptor 1
Background	similarity:Belongs to the LU7TM family.,