

GPR108 Polyclonal Antibody

Cat No: HR1AP3558

For research use only

Overview

Product Name	GPR108 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	GPR108
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR108; LUSTR2; Protein GPR108; Lung seven transmembrane receptor 2
Background	similarity:Belongs to the LU7TM family.,