

GPR113 Polyclonal Antibody

Cat No: HR1AP6093

For research use only

Overview

Product Name	GPR113 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	116kDa
GeneID?Human?	GPR113
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR113; PGR23; Probable G-protein coupled receptor 113; G-protein coupled receptor PGR23
Background	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily.,similarity:Contains 1 GPS domain.,