

## Frizzled-1 Polyclonal Antibody

Cat No: HR1AP3482

For research use only

### Overview

Product Name	Frizzled-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	71kDa
GeneID?Human?	FZD1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FZD1; Frizzled-1; Fz-1; hFz1; FzE1
Background	frizzled class receptor 1(FZD1) Homo sapiens Members of the 'frizzled' gene family encode 7-transmembrane domain proteins that are receptors for Wnt signaling proteins. The FZD1 protein contains a signal peptide, a cysteine-rich domain in the N-terminal extracellular region, 7 transmembrane domains, and a C-terminal PDZ domain-binding motif. The FZD1 transcript is expressed in various tissues. [provided by RefSeq, Jul 2008],