

CLPTM1 Polyclonal Antibody

Cat No: HR1AP5798

For research use only

Overview

Product Name	CLPTM1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	74kDa
GeneID?Human?	CLPTM1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CLPTM1; Cleft lip and palate transmembrane protein 1
Background	function:May play a role in T-cell development.,miscellaneous:A chromosomal translocation involving CLPTM1 is found in a family with cleft lip and palate. However, no etiologic link with the disease is characterized. Translocation t(2;19)(q11.2;q13.3).,similarity:Belongs to the CLPTM1 family.,tissue specificity:Widely expressed.,