

Frizzled-10 Polyclonal Antibody

Cat No: HR1AP5121

For research use only

Overview

Product Name	Frizzled-10 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	65kDa
GeneID?Human?	FZD10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FZD10; Frizzled-10; Fz-10; hFz10; FzE7; CD350
Background	frizzled class receptor 10(FZD10) Homo sapiens This gene is a member of the frizzled gene family. Members of this family encode 7-transmembrane domain proteins that are receptors for the Wingless type MMTV integration site family of signaling proteins. Most frizzled receptors are coupled to the beta-catenin canonical signaling pathway. Using array analysis, expression of this intronless gene is significantly up-regulated in two cases of primary colon cancer. [provided by RefSeq, Jul 2008],