

GPR152 Polyclonal Antibody

Cat No: HR1AP3566

For research use only

Overview

Product Name	GPR152 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	51kDa
GeneID?Human?	GPR152
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR152; PGR5; Probable G-protein coupled receptor 152; G-protein coupled receptor PGR5
Background	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,