

## T2R3 Polyclonal Antibody

Cat No: HR1AP7592

For research use only

### Overview

Product Name	T2R3 Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	TAS2R3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TAS2R3; Taste receptor type 2 member 3; T2R3
Background	<p>taste 2 receptor member 3(TAS2R3) Homo sapiens This gene encodes a member of a family of candidate taste receptors that are members of the G protein-coupled receptor superfamily and that are specifically expressed by taste receptor cells of the tongue and palate epithelia. These apparently intronless taste receptor genes encode a 7-transmembrane receptor protein, functioning as a bitter taste receptor. This gene is clustered with another 3 candidate taste receptor genes in chromosome 7 and is genetically linked to loci that influence bitter perception. [provided by RefSeq, Jul 2008],</p>