

## Smo Polyclonal Antibody

Cat No: HR1AP8354

For research use only

### Overview

Product Name	Smo Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	86kDa
GeneID?Human?	SMO
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SMO; SMOH; Smoothened homolog; SMO; Protein Gx
Background	smoothened, frizzled class receptor(SMO) Homo sapiens The protein encoded by this gene is a G protein-coupled receptor that interacts with the patched protein, a receptor for hedgehog proteins. The encoded protein transduces signals to other proteins after activation by a hedgehog protein/patched protein complex. [provided by RefSeq, Jul 2010],