

MC4-R Polyclonal Antibody

Cat No: HR1AP7314

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | MC4-R Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IHC-p,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 37kDa |
| GeneID?Human? | MC4R |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | MC4R; Melanocortin receptor 4; MC4-R |
| Background | melanocortin 4 receptor(MC4R) Homo sapiens The protein encoded by this gene is a membrane-bound receptor and member of the melanocortin receptor family. The encoded protein interacts with adrenocorticotrophic and MSH hormones and is mediated by G proteins. This is an intronless gene. Defects in this gene are a cause of autosomal dominant obesity. [provided by RefSeq, Jan 2010], |