

GPR22 Polyclonal Antibody

Cat No: HR1AP6709

For research use only

Overview

Product Name	GPR22 Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	GPR22
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR22; Probable G-protein coupled receptor 22
Background	G protein-coupled receptor 22(GPR22) Homo sapiens This gene is a member of the G-protein coupled receptor 1 family and encodes a multi-pass membrane protein. [provided by RefSeq, Jul 2008].