

GPR175 Polyclonal Antibody

Cat No: HR1AP5873

For research use only

Overview

Product Name	GPR175 Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	TPRA1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TPRA1; GPR175; PP6566; Transmembrane protein adipocyte-associated 1; Integral membrane protein GPR175
Background	TPRA1 (Transmembrane Protein Adipocyte Associated 1) is a Protein Coding gene. GO annotations related to this gene include G-protein coupled receptor activity.